DTSC Laboratory Update


Bruce LaBelle, June-Soo Park, Myrto Petreas Environmental Chemistry Laboratory DTSC

Report to Scientific Guidance Panel
Sacramento
November 2, 2010

Validated Methods


- Polybrominated Diphenyl Ethers (PBDEs)
- Perfluorinated Chemicals (PFCs)
- Polychlorinated Biphenyls (PCBs)
- Organochlorine Pesticides (DDT, DDE, HCB, HCH, Chlordane metabolites)
- Phenolic Compounds (e.g, OH-PCBs, OH-PBDEs, triclosan, pentachlorophenol)

CI-Pesticides J in CA Women





1960s (n=495); 1980s (n=30); 2008-09 (n=33)

PCBs ↓ **PBDEs** ↑ in **CA** Women



PCB/PBDE OH-Metabolites Are Potent Thyroid Endocrine Disruptors







PFC Analysis

-On-line automated SPE (Symbiosis) coupled with LC-MSMS (ABSciex 4000 Qtrap)

```
-QAQC:
Laboratory/Instrument Blanks 
MS control samples 
NIST SRM 
Inter-lab comparisons
CDC 
NYSL 
MNSL 
MNSL
```

PFCs in California Women





Green bar: NHANES (1999-00) and (2003-04)

CDC Cooperative Agreement

- Personnel:
 - Dr. Suhash Harwani
 - Dr. Tan Guo
 - -Third Candidate in Process

- Collaboration with CDPH on QA/QC, Sample Management, LIMS
- **Equipment: LC-MSMS**

New/Alternative BFRs: Soon to come...

	GC-MS Method	Detect in Human Serum?	Further method Validation?
ATE	OK	n.d.?	Yes (GC-HRMS)
a,b-TBECH	OK	n.d.?	Yes (GC-HRMS)
PBT	OK	n.d.?	Yes (GC-HRMS)
PBEB	OK	n.d.?	Yes (GC-HRMS)
DPTE	OK	n.d.?	Yes (GC-HRMS)
HBB	OK	detect	Yes (GC-HRMS)
DBDPE	OK	n.d.?	Yes (GC-HRMS)

Soon to Come...

Develop off-line SPE/LC-MSMS methods for phenolic compounds (e.g. PCB/PBDE OH-metabolites, triclosan, BPA, etc.)

Analyze serum samples from California study populations