

Report to Scientific Guidance Panel:

CDPH

Environmental Health Laboratory Update

Jianwen She, PhD

Environmental Health Laboratory
CDPH

Oakland, May 24, 2010

Update Overview

- New staff
- Laboratory set up and instrumentation
- New method development
- Method performance
- Proposed year 2 activities

New Cooperative Agreement Staff

- Dr. Binfeng Xia, Laboratory Scientist
- Dr. Dongli Wang, Laboratory Scientist
- One Laboratory Information Management Specialist to be hired in July 2010

These are in addition to 4 lab staff hired through February 2010 (2 ELS plus QMC and SMS)

Laboratory Set up

- Currently in use:
 - Sample receiving and handling room (D-366)
 - Organics and inorganics analytical lab (D-306)
- Identified for future use:
 - Organics analytical lab (D-263 est. August 2010)
 - Organics analytical lab (D-262 est. Dec 2010)



New Instruments

- Instruments ordered:
 - LC-MS/MS (AB QTRAP 5500): OP specific metabolites and pyrethroid metabolites.
 - LC-MS/MS (Agilent 6460): environmental phenols
 - GC-MS/MS (Agilent 7000): OP common metabolites: dialkyl phosphate (DAPs)
- Still to be purchased:
 - IC-MS/MS: perchlorate
 - Sample preparation equipment: lab automation



Method Development

Validated	In Progress	Exploratory
1. Blood metals (Pb, Cd, Hg, Mn)	1. Metals in urine	1. DAPs
2. Chlorpyrifos (TCPy) Pyrethroid (3-PBA)	2. Arsenic speciation	2. Perchlorate
3. Phthalate (MbP, MeP)	3. Bisphenol – A	
4. PAH (3-Phen)		
5. Creatinine		

Method Precision: 3-Phen



Method Accuracy: Blood Metals



Year 2 Activities

- **Expand upon existing methods:**
 - Hydroxy-PAH (currently 1 analyte; will increase in Year 2)
 - Phthalate metabolites (currently 2 analytes; will increase in Year 2)
 - OP specific metabolite (currently 1 analyte; will increase in Year 2)
- **Methods in progress:**
 - Metals in urine (Oct. 2010)
 - Arsenic speciation (June 2011)
 - Bisphenol A (June 2011)
- **Increase capacity**
 - Procedure automation
 - Enhance through-put

